



Myeloperoxidase (MPO) ELISA Kit Instructions

For the quantitative determination of MPO in
human serum and plasma

**Catalog #80727
96 Assays**

For research use only. Not for use in diagnostic procedures.

Crystal Chem, Inc.
955 Busse Road
Elk Grove Village, IL 60007, USA
Tel: (630) 889-9003 Fax: (630) 889-9021
E-mail: sales@crystalchem.com
Order online: www.crystalchem.com

TABLE OF CONTENTS

| | |
|--|---|
| A. <i>Intended Use</i> | 1 |
| B. <i>Introduction</i> | 1 |
| C. <i>Principles of the Assay</i> | 1 |
| D. <i>Kit Storage</i> | 1 |
| E. <i>Assay Materials</i> | |
| E.1. Materials provided | 1 |
| E.2. Materials required but not provided | 1 |
| F. <i>Assay Precautions</i> | 2 |
| G. <i>Maximizing Kit Performance</i> | 2 |
| H. <i>Sample Collection</i> | 2 |
| I. <i>Assay Procedure</i> | |
| I.1. Preparation of reagents | 2 |
| I.2. Preparation of working standards | 3 |
| I.3. Dilution of samples..... | 3 |
| I.4. Assay procedure | 4 |
| I.5. Determining the MPO concentration..... | 4 |
| J. <i>Performance characteristics</i> | |
| J.1. Assay range | 4 |
| J.2. Precision | 4 |
| Warranty..... | 5 |

A. Intended Use

The Myeloperoxidase (MPO) ELISA kit is for the quantitative determination of MPO in human serum and plasma. Please read the complete kit insert before performing this assay. The kit is for RESEARCH USE ONLY. It is not intended for use in diagnostic procedures. For in Vitro Use Only.

B. Introduction

Myeloperoxidase (MPO) is a peroxidase enzyme synthesized in neutrophil granulocytes and monocyte precursor cells. It is a 150kDa protein comprised of two light chains and two heavy chains. Recent reports suggest that MPO levels might serve as predictors for myocardial infarctions or markers for coronary artery disease.

C. Principle of the Assay

The Myeloperoxidase (MPO) ELISA kit is a double antibody sandwich ELISA. In the assay, the MPO present in the sample binds with anti-MPO antibodies adsorbed to the surface of the microplate. In the following step after washing, a specific anti-MPO-Antibody binds, in turn, to the immobilized MPO. The second antibody is conjugated to HRP, and in the closing substrate reaction, the MPO levels of the samples can be measured by color intensity.

D. Kit Storage

1. Upon receipt of the Myeloperoxidase (MPO) ELISA kit, store it at 2-8°C and avoid light exposure (do not freeze the kit or hold it at temperatures above 25°C).
2. The kit should not be used after the expiration date.

E. Assay Materials**E.1. Materials provided****TABLE 1 CONTENTS OF KIT**

| Mark | Description | Amount |
|---------|---------------------------------------|---------------|
| MIC | Antibody-coated Microplate (12 x 8) | 1 pack |
| CAL | MPO Calibrator (Lyophilized) | 1 x 1 vial |
| DIL | Diluent (5X Concentrate) | 1 x 50 mL |
| AB CONJ | Antibody Conjugate (100X Concentrate) | 1 vial/150 µL |
| WASH | Wash Buffer (20X Concentrate) | 1 x 50 mL |
| SUB | Substrate Solution | 1 x 12 mL |
| STOP | Stop Solution | 1 x 12 mL |

E.2. Materials required but not provided

Micropipettes and disposable tips
 Standard laboratory glassware for buffer and reagent preparation
 Distilled or deionized water
 Polypropylene microtubes
 Microplate reader (capable of reading A_{450} and A_{630} values)

F. Assay Precautions

1. Only appropriately-trained personnel should use the kit. Laboratory personnel should wear suitable protective clothing. All chemicals and reagents should be considered potentially hazardous. Avoid ingestion and contact with skin and eyes. In case of contact with eyes or skin, flush immediately with water and contact a medical professional.
2. Some assay components may contain human sourced materials. Accordingly, all assay components should be handled as if potentially infectious using safe laboratory procedures.
3. Do not use the reagents after the expiration date.
4. Reagents are light sensitive and should be protected from sunlight.

G. Maximizing Kit Performance

1. Each standard and sample should be assayed in duplicate.
2. Given the sample volumes required (30 μ L), pipetting should be done as carefully as possible. A high quality 50 μ L or better precision pipette should be used for such volumes. Drops of liquid adhering to the outside of the pipette tips should be removed by wiping to ensure the highest degree of accuracy.
3. In order to prevent the microplate wells from drying out and to get the best results, samples and reagents should be dispensed quickly into the wells.
4. The same sequence of pipetting and other operations should be maintained in all procedures.
5. Do not mix reagents that have different lot numbers.

H. Sample Collection

Blood should be collected by venipuncture. The serum should be separated from the cells after clot formation by centrifugation. For plasma samples, blood should be collected into a container with an anticoagulant and then centrifuged. The samples should be assayed immediately or aliquoted and stored at -20°C. Avoid repeated freeze-thaw cycles. Samples with excessive hemolysis should not be used. *Note: Azide and thimerosal at concentrations higher than 0.1% inhibit the enzyme reaction.*

I. Assay Procedure

I.1. Preparation of reagents

1. Antibody-coated microplate
Provided as ready to use. Protect from moisture.
2. Calibrator
The calibrator needs to be reconstituted with 1.0 mL of distilled or deionized water to prepare the 0.92 μ g/mL calibration standard. The reconstituted standard should be immediately used, and the remaining portion should be aliquoted and stored frozen. Avoid repeated freeze/thaw cycles. The reconstituted standard is stable for 1 hour, and it should be diluted with 1X Diluent to prepare working standards as described in I.2.
3. Diluent (5X Concentrated)
The diluent has to be diluted 1:5 with distilled or deionized water prior to use. For example, 50 mL of diluent must be diluted with 200 mL of distilled or deionized water. Diluent is stable for at least one week after dilution stored at 2-8°C. Prepare only as needed.

4. Antibody Conjugate (100X Concentrated)

The Antibody Conjugate has to be diluted 1:100 with 1X Diluent prior to use. For example, 10 μ L of the Antibody Conjugate must be diluted with 990 μ L of 1X Diluent. The working Antibody Conjugate is stable for up to 1 hour when stored in the dark. Accordingly, Antibody Conjugate should be prepared only as needed just prior to use.

5. Wash Buffer (20X Concentrated)

The wash buffer has to be diluted 1:20 with distilled or deionized water prior to use. For example, 50 mL of wash buffer must be diluted with 950 mL of distilled or deionized water. Wash buffer is stable for at least one week after dilution at room temperature or at 2-8°C.

6. Substrate Solution

Provided as ready to use.

7. Stop Solution

Provided as ready to use.

I.2. Preparation of working standards

1. Pipette 800 μ L of 1X diluent and 120 μ L of the standard (0.92 μ g/mL) into a polypropylene microtube labeled 120 ng/mL, and mix thoroughly.
2. Pipette 300 μ L of 1X diluent into six polypropylene microtubes labeled 60, 30, 15, 7.50, 3.75, and 1.88 ng/mL.
3. Dispense 300 μ L of the 120 ng/mL standard into the 60 ng/mL microtube, and mix thoroughly.
4. Dispense 300 μ L of the 60 ng/mL standard into the 30 ng/mL microtube, and mix thoroughly.
5. Dispense 300 μ L of the 30 ng/mL standard into the 15 ng/mL microtube, and mix thoroughly.
6. Dispense 300 μ L of the 15 ng/mL standard into the 7.50 ng/mL microtube, and mix thoroughly.
7. Dispense 300 μ L of the 7.50 ng/mL standard into the 3.75 ng/mL microtube, and mix thoroughly.
8. Dispense 300 μ L of the 3.75 ng/mL standard into the 1.88 ng/mL microtube, and mix thoroughly.
9. Dispense 600 μ L of 1X diluent into one polypropylene microtube labeled 0 ng/mL. You now have working standards of 120, 60, 30, 15, 7.50, 3.75, 1.88, and 0 ng/mL.

Please note: Working standards should be prepared immediately prior to use. Avoid foaming.

I.3. Dilution of samples

Samples need to be diluted with 1X diluent for use with the assay. A sample dilution of 1:10 using 30 μ L of sample is generally suitable. Dilute the samples immediately before use, and use within 1 hour after dilution.

Serum/Plasma Samples: To prepare the 1:10 dilution, mix 30 μ L of sample with 270 μ L of 1X Diluent.

Since MPO levels can vary, dilution ratios may need to be adjusted as appropriate.

I.4. Assay procedure

Prior to running the assay, all reagents should be brought to room temperature for at least 30 minutes. Reagents should be stored at 2-8°C immediately after use. Before use, mix the reagents thoroughly by gentle agitation or swirling.

1. In each well, add 100 μ L of diluted sample or working standards.
2. Incubate plate for 60 mins at room temperature. Keep plate covered and level.
3. Aspirate well contents and wash four times using 300 μ L of 1X Wash Buffer per well. After each wash, remove any remaining solution by inverting and tapping the plate firmly on a clean paper towel.
4. In each well, add 100 μ L of working Antibody Conjugate.
5. Incubate plate for 20 mins at room temperature. Keep plate covered in the dark and level.
6. Aspirate well contents and wash four times using 300 μ L of 1X Wash Buffer per well. After each wash, remove any remaining solution by inverting and tapping the plate firmly on a clean paper towel.
7. Add 100 μ L of the Substrate Solution in each well.
8. Incubate plate for 10 mins in dark room at room temperature.
9. Stop the reaction by adding 100 μ L of Stop Solution.
10. Measure absorbance as soon as possible, and within 30 mins, using a plate reader (measure A_{450} values and subtract A_{650} values).

I.5. Determining the MPO concentration

1. Using computer software, construct the MPO calibration curve by plotting the mean absorbance value for each standard (not blank) on the Y axis versus the corresponding MPO concentration on the X axis. A four parametric logistic (4-PL) curve fit or second order polynomial (quadratic) are suitable for the evaluation.

Note: A calibration curve should be plotted every time the assay is performed.

2. MPO concentrations in the samples are interpolated using the calibration curve and mean absorbance values for each sample. For diluted samples, the values obtained must be multiplied by the dilution factor to obtain the final MPO concentration (expressed in ng/mL).

Note: Samples with high MPO concentrations (ie. fall above the range of the assay) should be further diluted and rerun.

J. Performance characteristics

J.1. Assay range

The Myeloperoxidase (MPO) ELISA Kit has an assay range from 1.88 - 120 ng/mL.

J.2. Precision

The assay has a within-run and total precision of CV < 20%.

Warranty

Crystal Chem Inc. makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Buyer assumes all risk and liability resulting from the use of this product.

THERE IS NO WARRANTY OF MERCHANTABILITY OF THE PRODUCTS, OR THAT SUCH PRODUCTS ARE FIT FOR ANY PARTICULAR PURPOSE. CRYSTAL CHEM INC'S LIABILITY SHALL NOT EXCEED THE RETURN OF THE PURCHASE PRICE, AND UNDER NO CIRCUMSTANCES SHALL CRYSTAL CHEM INC. BE LIABLE FOR SPECIAL OR CONSEQUENTIAL DAMAGES, OR EXPENSES ARISING DIRECTLY OR INDIRECTLY FROM THE USE OF THIS PRODUCT.

Myeloperoxidase (MPO) ELISA
Kit Instructions

For the quantitative determination of MPO in
human serum and plasma

Catalog #80727
25 assays