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1. Riggs L. Causes of age-related bone loss and fractures, Osteoporosis: Physiological Basis, Assessment, and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 7-16, 1990.
2. Status of Research in Osteoporosis. Osteoporosis Research, Education and Health Promotion. U.S. Department of Health and Human Services, Public Health Service, National Institutes of Health.
3. Canalis E. Regulation of bone remodeling, Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. Favus MJ, ed., Raven Press (New York), pp. 33-37, 1993.
4. Kleerekoper M and Avioli LV. Evaluation and treatment of postmenopausal osteoporosis, Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. Favus MJ, ed., Raven Press (New York), pp. 223-229, 1993.
5. Lindsay R. Osteoporosis. A Guide to Diagnosis, Prevention, and Treatment. Raven Press (New York), pp. 1-39, 1992.
6. Epstein S. Serum and urinary markers of bone remodeling: assessment of bone turnover. Endocrine Reviews 9:437-448, 1988.
7. Whyte MP. Alkaline phosphatase and the measurement of bone formation, Clinical Disorders of Bone and Mineral Metabolism. Potts JT and Frame B, eds., pp. 120-125, 1983.
8. Garnero P and Delmas P. Assessment of the serum levels of bone alkaline phosphatase with a new immunoradiometric assay in patients with metabolic bone disease. J Clin Endocrin Metab 77:1046-1053, 1993.
9. Cooper EH, Whelan P, and Purves D. Bone alkaline phosphatase and prostate-specific antigen in the monitoring of prostate cancer. The Prostate 25:236-242, 1994.
10. Van Straalen JP, Sanders E, Prummel MF, and Sanders GTB. Bone alkaline phosphatase as indicator of bone formation. Clin Chim Acta 201:27-33, 1991.
11. Merkow RL and Lane JM. Metabolic bone disease and Paget's disease in the elderly: II. Paget's Disease. Clin Rheum Dis 12:70, 1986.
12. Lindsay R. Sex steroids in the pathogenesis and prevention of osteoporosis, Osteoporosis: Etiology, Diagnosis, and Management. Riggs L and Melton LJ, eds., Raven Press (New York), pp. 333-358, 1988.
13. Peck W. The pathogenesis of postmenopausal osteoporosis, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 3-6, 1990.
14. Christiansen C. Osteoporosis: Diagnosis and management today. Bone 17:513S-516S, 1995.
15. Turner R, Riggs B and Spelsberg T. Skeletal effects of estrogen. Endocrine Reviews 15:275-300, 1994.
16. Kanis J. Pathogenesis of osteoporosis and fracture, Osteoporosis. London: Blackwell Science Ltd.; pp. 22-55, 1994.
17. Mack TM and Ross RK. A current perception of HRT risks and benefits, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 161-178, 1990.
18. Riggs L and Melton LJ. The prevention and treatment of osteoporosis. N Engl J Med 327:620-627,

1992.

19. Guidelines for Preclinical and Clinical Evaluation of Agents Used in the Prevention or Treatment of Postmenopausal Osteoporosis. Division of Metabolism and Endocrine Drug Products, Food and Drug Administration. April, 1994.
20. Ribot C, Tremollieres F and Pouilles JM. Prevention of postmenopausal bone loss by long-term parenteral administration of 17β estradiol: comparison of percutaneous and transdermal route, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 137-145, 1990.
21. Delmas P. Biochemical markers of bone turnover. Methodology and clinical use in osteoporosis. Am J Med 91:59S-63S, 1991.
22. Duda RJ, O'Brien JF, Katzman JA, Peterson JM, Mann KG, and Riggs BL. Concurrent assays of circulating bone gla-protein and bone alkaline phosphatase: effects of sex, age, and metabolic bone disease. J Clin Endocrin Metab 66:951-957, 1988.
23. HHS Publication No. (CDC) 93-8395. 3rd edition. May 1993. Biosafety in Microbiological and Biomedical Laboratories. Washington, DC: US Government Printing Office, 1993.
24. Manual Guide - Safety Management No. CDC-22. Decontamination of Laboratory Sink Drains to Remove Azide Salts. April 30, 1976. Atlanta, GA: Centers for Disease Control.
25. Hill CS and Wolfert RL. The preparation of monoclonal antibodies which react preferentially with human bone alkaline phosphatase and not liver alkaline phosphatase. Clin Chim Acta 186:315-320, 1989.
26. Garnero P, Shih WJ, Gineyts E, Karpf DB, and Delmas PD. Comparison of new biochemical markers of bone turnover in late postmenopausal osteoporotic women in response to alendronate treatment. J Clin Endocrin Metab 79:1693-1700, 1994.
27. Liberman UA, Weiss SR, Broll J, Minne HW, Quan H, Bell NH, Rodriguez-Portales J, Downs RW, Dequeker J, Favus M, Seeman E, Recker RR, Capizzi T, Santora AC, Lombardi A, Shah RV, Hirsch LJ and Karpf DB. For the Alendronate Phase III Osteoporosis Treatment Study Group. Effect of oral alendronate on bone mineral density and the incidence of fractures in postmenopausal osteoporosis. N Engl J Med 333:1437-1443, 1995.
28. Dawson-Hughes B, Dallal GE, Krall EA, Sadowski L, Sahyoun N and Tannenbaum S. A controlled trial of the effect of calcium supplementation on bone density in postmenopausal women. N Engl J Med 323:878-883, 1990.
29. Reid IR, Ames RW, Evans MC, Gamble GD, and Sharpe SJ. Long-term effects of calcium supplementation on bone loss and fractures in postmenopausal women: a randomized controlled trial. Am J Med 98:331-335, 1995.
30. Devogelaer JP, Broll H, Correa-Rotter R, Cumming DC, De Deuchchaisnes CN, Geusens P, Hosking D, Jaeger P, Kaufman JM, Leite M, Leon J, Liberman U, Menkes CJ, Meunier PJ, Reid I, Rodriguez J, Romanowicz A, Seeman E, Vermeulen A, Hirsch LJ, Lombardi A, Plezia K, Santora AC, Yates AJ and Yuan W. Oral alendronate induces progressive increases in bone mass of the spine, hip, and total body over 3 years in postmenopausal women with osteoporosis. Bone 18:141-150, 1996.
31. Overgaard K, Alexandersen P, Riis BJ and Christiansen C. Evaluation of a new commercial IRMA for bone-specific alkaline phosphatase during treatment with hormone replacement therapy and calcitonin. Clin Chem 42:973-974, 1996.
32. Soletormos G, Schioler V, Nielsen D, Skovgaard T, and Dombrowsky P. Interpretation of results for tumor markers on the basis of analytical imprecision and biological variation. Clin Chem 10:2077-2083, 1993.
33. Blumsohn A, Hannon RA, Al-Dehaimi AW, and Eastell R. Short-term intraindividual variability of markers of bone turnover in healthy adults. J Bone Min Res 9:1533, 1994.
34. Raisz L, Wiita B, Artis A, Bowen A, Bowen A, Schwartz S, Trahiotis M, Shoukri K and Smith J. Comparison of the effects of estrogen alone and estrogen plus androgen on biochemical markers of bone formation and resorption in postmenopausal women. J Clin Endocrin Metab 81:37-43, 1996.
35. Kanis J. Treatment of generalized osteoporosis with inhibitors of bone resorption, Osteoporosis. London: Blackwell Science Ltd.; pp. 168-195, 1994.
36. Kennedy J, Carey R, Coolen R, Garber C, Hartman A, Lee H, Leitz V, Levine J, McLean M, Osberg I, Steindel S, and Sylvester E. User Evaluation of Precision Performance of Clinical Chemistry Devices. Tentative Guideline, NCCLS 4:8 EP5-T2, 1984.
37. Moss D and Whitby L. A simplified heat-inactivation method for investigating alkaline phosphatase isoenzymes in serum. Clin Chem Acta 61:63, 1975.
38. Price C, Mitchell C, Moriarty J, Gray M and Noonan K. Mass versus activity of an immunometric assay for bone alkaline phosphatase in serum. Ann Clin Biochem 32:405, 1995.

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注册人/生产企业名称: Immunodiagnostic Systems Limited (IDS Ltd)

住所: 10 Didcot Way, Boldon Business Park, Boldon, Tyne & Wear, NE35 9PD, UK

生产地址: 10 Didcot Way, Boldon Business Park, Boldon, Tyne & Wear, NE35 9PD, UK

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网址: www.idsplc.com